

Product Datasheet

Tricyclic Antidepressant (TCA) Antibody PP141IF

Product Name	Tricyclic Antidepressant (TCA) Antibody
Catalog Number	PP141IF
Species/Host	Sheep
Tested Applications	ELISA
Immunogen	TCA-BTG
Brand	EastCoast Bio
Form/Appearance	DEAE Purified Polyclonal Antibody, Ig Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Preservatives	0.05% Sodium Azide