

## Product Datasheet

### Tricyclic Antidepressant (TCA) Antibody PP142IF

<b>Product Name</b>	Tricyclic Antidepressant (TCA) Antibody
<b>Catalog Number</b>	PP142IF
<b>Species/Host</b>	Sheep
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	TCA-BgG
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	DEAE Purified Polyclonal Antibody, Ig Fraction
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.0
<b>Preservatives</b>	0.05% Sodium Azide