

Product Datasheet

Cotinine-3 Antibody PP192IF

Product Name Cotinine-3 Antibody

Catalog Number PP192IF

Species/Host Rabbit

Tested Applications ELISA

Immunogen Cotinine-3-BgG

Brand EastCoast Bio

Form/Appearance Protein A Purified Polyclonal Antibody, Ig Fraction

Storage Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Clonality Polyclonal

Buffer 10 mM Phosphate Buffered Saline, pH 7.0

Matched Pairs P82-99-01A (capture) or P92-99-01H (detection)

Preservatives 0.05% Sodium Azide