

Product Datasheet

Cotinine-3 Antibody PP192IF

Product Name	Cotinine-3 Antibody
Catalog Number	PP192IF
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	Cotinine-3-BgG
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Polyclonal Antibody, Ig Fraction
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Matched Pairs	P82-99-01A (capture) or P92-99-01H (detection)
Preservatives	0.05% Sodium Azide