

## Product Datasheet

### Cotinine-3 Antibody PP193IF

<b>Product Name</b>	Cotinine-3 Antibody
<b>Catalog Number</b>	PP193IF
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Cotinine-3-BgG
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Polyclonal Antibody, Ig Fraction
<b>Storage</b>	Short Term ( $\leq$ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.0
<b>Matched Pairs</b>	P82-99-01A (capture) or P92-99-01H (detection)
<b>Preservatives</b>	0.05% Sodium Azide