

Product Datasheet

Cotinine-4 Antibody PP198IF

Product Name	Cotinine-4 Antibody
Catalog Number	PP198IF
Species/Host	Sheep
Tested Applications	LF, ELISA
Immunogen	Cotinine-4-BgG
Brand	EastCoast Bio
Form/Appearance	DEAE Purified Polyclonal Antibody, Ig Fraction
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>90%
Buffer	20 mM Sodium Phosphate, pH 7.1 with 150 mM Sodium Chloride
Matched Pairs	P82-99-02A (capture) or P92-99-02H (detection)