

## Product Datasheet

### Ketamine Antibody PP211IF

**Product Name** Ketamine Antibody

**Catalog Number** PP211IF

**Species/Host** Rabbit

**Tested Applications** ELISA

**Immunogen** Ketamine - BgG

**Brand** EastCoast Bio

**Form/Appearance** Protein A Purified Polyclonal Antibody, Ig Fraction

**Storage** Short Term ( $\leq 2$  weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Clonality** Polyclonal

**Buffer** 10 mM Phosphate Buffered Saline, pH 7.0

**Matched Pairs** P82-99-31A (capture) or P92-99-31H (detection)

**Preservatives** 0.05% Sodium Azide