

Product Datasheet

Varicella-Zoster Virus (VZV-gE) Antibody VA315

Product Name	Varicella-Zoster Virus (VZV-gE) Antibody
Catalog Number	VA315
Species/Host	Mouse
Tested Applications	ELISA, WB, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	gE