

Product Datasheet

Thyroxine (T4) Antibody Y306

Product Name	Thyroxine (T4) Antibody
Catalog Number	Y306
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	T4-KLH
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (\leq 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate, pH 7.4 with 150 mM NaCl
Specificity	T4
Preservatives	0.05% Sodium Azide