

Product Datasheet

Recombinant Human CD25 - IL-2RA (CHO) 705-H27-100

Product Name	Recombinant Human CD25 - IL-2RA (CHO)
Catalog Number	705-H27-100
Description	Recombinant Protein
Species/Host	CHO Cells
Alternative Names	Interleukin-2 Receptor subunit Alpha
Brand	Diaclone
Form/Appearance	Liquid ; Phosphate-buffered saline pH7.4. Sterile-filtered through 0.22 µm.
Storage	Store at minimum -20°C. Avoid repeated freeze-thaw cycles.
Purity	>90% by SDS PAGE
Sequence	Glu22-Gln240 (ECD)
Tag	HIS Tag C-Terminus