

Product Datasheet

Recombinant Human CD25 - IL-2RA (CHO) 705-H27-100

| | |
|--------------------------|-----------------------------------------------------------------------------|
| Product Name | Recombinant Human CD25 - IL-2RA (CHO) |
| Catalog Number | 705-H27-100 |
| Description | Recombinant Protein |
| Species/Host | CHO Cells |
| Alternative Names | Interleukin-2 Receptor subunit Alpha |
| Brand | Diaclone |
| Form/Appearance | Liquid ; Phosphate-buffered saline pH7.4. Sterile-filtered through 0.22 µm. |
| Storage | Store at minimum -20°C. Avoid repeated freeze-thaw cycles. |
| Purity | >90% by SDS PAGE |
| Sequence | Glu22-Gln240 (ECD) |
| Tag | HIS Tag C-Terminus |