

## Product Datasheet

### Recombinant Human SARS-CoV-2 Spike RBD Variant Mu B.1.621 (HEK) 715-H26-0BU

<b>Product Name</b>	Recombinant Human SARS-CoV-2 Spike RBD Variant Mu B.1.621 (HEK)
<b>Catalog Number</b>	715-H26-0BU
<b>Description</b>	Recombinant Protein
<b>Species/Host</b>	HEK cells
<b>Alternative Names</b>	Variant Mu, Colombia, 21H
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Liquid ; Phosphate-buffered saline pH7.4. Sterile-filtered through 0.22 µm.
<b>Storage</b>	Store at minimum -20°C. Avoid repeated freeze-thaw cycles.
<b>Purity</b>	>95% by SDS PAGE
<b>Biological Activity</b>	Recognized by CR3022 anti RBD recombinant antibody
<b>Sequence</b>	Arg319-Phe541 with R346K; E484K; N501Y mutations
<b>Tag</b>	His Tag C-Terminus