

Product Datasheet

Anti-Human CD226 Monoclonal Antibody, Azide Free Clone B-G51 857.360

Product Name Anti-Human CD226 Monoclonal Antibody, Azide Free Clone B-G51

Catalog Number 857.360

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, FC

Immunogen Recombinant human CD226

Alternative Names DNAM-1

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile filtered through 0.22µm. Carrier and preservative free.

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG2a Kappa light chain

Clonality Monoclonal

Specificity Recognizes both natural and recombinant human CD226, a 65 kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells