

Product Datasheet

Anti-Human CD226 Monoclonal Antibody, Unconjugated Clone B-G51 857.363.020

Product Name	Anti-Human CD226 Monoclonal Antibody, Unconjugated Clone B-G51
Catalog Number	857.363.020
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	Recombinant human CD226
Alternative Names	DNAM-1
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognizes both natural and recombinant human CD226, a 65 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells