

Product Datasheet

2'-O-Me sensitive RT-KTQ DNA polymerase mutant 9999

Product Name	2'-O-Me sensitive RT-KTQ DNA polymerase mutant
Catalog Number	#9999
Description	5 μ M, 20 μ l For more information, please check : 2'-O-Me sensitive RT-KTQ DNA polymerase mutant
Tested Applications	End-Point, Real-Time
Brand	myPOLS Biotec
Storage	-20°C
References	Aschenbrenner, J., & Marx, A. (2016). Direct and site-specific quantification of RNA 2'-O-methylation by PCR with an engineered DNA polymerase. Nucleic acids research, 44(8), 3495-3502.
Content	S pack: 5 μ M, 1 x 20 μ l 2'-O-Me sensitive RT-KTQ DNA polymerase mutant M pack: 5 μ M, 1 x 100 μ l 2'-O-Me sensitive RT-KTQ DNA polymerase mutant