

## Product Datasheet

### RT-KTQ2 DNA polymerase 9998

<b>Product Name</b>	RT-KTQ2 DNA polymerase
<b>Catalog Number</b>	#9998
<b>Description</b>	5 $\mu$ M, 20 $\mu$ l For more information, please check : <b>RT-KTQ2 DNA polymerase</b>
<b>Tested Applications</b>	End-Point, Real-Time
<b>Brand</b>	myPOLS Biotec
<b>Storage</b>	-20°C
<b>References</b>	Blatter, N., Bergen, K., Nolte, O., Welte, W., Diederichs, K., Mayer, J., ... & Marx, A. (2013). Structure and Function of an RNA-Reading Thermostable DNA Polymerase. <i>Angewandte Chemie International Edition</i> , 52(45), 11935-11939.
<b>Content</b>	S pack: 5 $\mu$ M, 1 x 20 $\mu$ l RT-KTQ2 DNA polymerase M pack: 5 $\mu$ M, 1 x 100 $\mu$ l RT-KTQ2 DNA polymerase