

## Product Datasheet

### m6A sensitive RT-KTQ DNA polymerase mutant 9997

<b>Product Name</b>	m6A sensitive RT-KTQ DNA polymerase mutant
<b>Catalog Number</b>	#9997
<b>Description</b>	5 $\mu$ M, 20 $\mu$ l For more information, please check : <b>m6A sensitive RT-KTQ DNA polymerase mutant</b>
<b>Tested Applications</b>	End-Point, Real-Time
<b>Brand</b>	myPOLS Biotec
<b>Storage</b>	-20°C
<b>References</b>	Aschenbrenner, J., Werner, S., Marchand, V., Adam, M., Motorin, Y., Helm, M., & Marx, A. (2018). Engineering of a DNA polymerase for direct m6A sequencing. <i>Angewandte Chemie International Edition</i> , 57(2), 417-421.
<b>Content</b>	S pack: 5 $\mu$ M, 1 x 20 $\mu$ l m6A sensitive RT-KTQ DNA polymerase mutant  M pack: 5 $\mu$ M, 1 x 100 $\mu$ l m6A sensitive RT-KTQ DNA polymerase mutant