

Product Datasheet

m6A sensitive RT-KTQ DNA polymerase mutant 9997

| | |
|----------------------------|--|
| Product Name | m6A sensitive RT-KTQ DNA polymerase mutant |
| Catalog Number | #9997 |
| Description | 5 μ M, 20 μ l For more information, please check : m6A sensitive RT-KTQ DNA polymerase mutant |
| Tested Applications | End-Point, Real-Time |
| Brand | myPOLS Biotec |
| Storage | -20°C |
| References | Aschenbrenner, J., Werner, S., Marchand, V., Adam, M., Motorin, Y., Helm, M., & Marx, A. (2018). Engineering of a DNA polymerase for direct m6A sequencing. <i>Angewandte Chemie International Edition</i> , 57(2), 417-421. |
| Content | S pack: 5 μ M, 1 x 20 μ l m6A sensitive RT-KTQ DNA polymerase mutant M pack: 5 μ M, 1 x 100 μ l m6A sensitive RT-KTQ DNA polymerase mutant |