

## Product Datasheet

### Cotinine-4-HRP Antigen 170049

<b>Product Name</b>	Cotinine-4-HRP Antigen
<b>Catalog Number</b>	170049
<b>Description</b>	Synthetic Cotinine, HRP conjugated
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid; Antigen-HRP conjugate
<b>Storage</b>	Storage at 2-8 °C. Do not freeze.
<b>Note</b>	This product was previously available with a Catalog code P92-99-02H. It has been produced using the same conjugation protocol.
<b>Buffer</b>	1:1 0.1 M Tris base / FBS, pH 7.0, 0.0016% Gentamicin sulfate, 0.005% thimerosal
<b>Matched Pairs</b>	PP194IF, PP195IF, PP198IF (capture)
<b>Shelf Life</b>	Unspecified