

Product Datasheet

Benzoylecgonine (BE) Antibody HM1207

Product Name	Benzoylecgonine (BE) Antibody
Catalog Number	HM1207
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Benzoylecgonine and Cocaine
Matched Pairs	LA004 (capture)
Preservatives	0.05% Sodium Azide