

Product Datasheet

Progesterone 17-OH Antibody HM1226

| | |
|----------------------------|--|
| Product Name | Progesterone 17-OH Antibody |
| Catalog Number | HM1226 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | 17-OH Progesterone-KLH |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG2b |
| Clonality | Monoclonal |
| Purity | >90% by SDS-PAGE |
| Buffer | 0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride |
| Cross Reactivity | 17-OH Hydroxyprogesterone >100%; Progesterone, Estrone, 17 α -Ethinylestradiol, 17 β -Estradiol, Estriol, Testosterone, Testosterone glucuronide, 5 α -Dihydrotestosterone, Epitestosterone, 6 β -Hydroxytestosterone, 17 α -Methyltestosterone, Aldosterone, Androstenedione <1% |
| Specificity | 17-OH Hydroxyprogesterone |
| Matched Pairs | LA330 (capture) |
| Preservatives | 0.05% Sodium Azide |