

Product Datasheet

Hantaan Virus (NP) Antibody HM1241

| | |
|----------------------------|---|
| Product Name | Hantaan Virus (NP) Antibody |
| Catalog Number | HM1241 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | Recombinant Hantaan Virus Nucleoprotein |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Purity | >90% by SDS-PAGE |
| Buffer | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride |
| Preservatives | 0.05% Sodium Azide |