

Product Datasheet

Zika Virus (ZIKV-NS1) Antibody HM1216

Product Name	Zika Virus (ZIKV-NS1) Antibody
Catalog Number	HM1216
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Native Zika Virus
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with Dengue, Chikungunya, Yellow Fever, Japanese Encephalitis, Tick-borne Encephalitis, or West Nile Virus
Specificity	Zika Virus NS1
Matched Pairs	HM1215 or HM334 (detection)
Preservatives	0.05% Sodium Azide