

Product Datasheet

Digoxin Antibody HM1210

Product Name	Digoxin Antibody
Catalog Number	HM1210
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Digoxin-BSA
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl
Specificity	Digoxin
Preservatives	0.05% Sodium Azide