

Product Datasheet

6-FAM (6-Carboxyfluorescein) Antibody HM1233

Product Name	6-FAM (6-Carboxyfluorescein) Antibody
Catalog Number	HM1233
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	6-FAM-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	0.05M Borate, pH 9.0
Specificity	6-Carboxyfluorescein
Preservatives	0.05% Sodium Azide