

Product Datasheet

FITC (Fluorescein isothiocyanate) Antibody HM1239

Product Name	FITC (Fluorescein isothiocyanate) Antibody
Catalog Number	HM1239
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	FITC-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Specificity	FITC
Preservatives	0.05% Sodium Azide