

Product Datasheet

Neurofilament H (NfH) Antibody HM1247

Product Name	Neurofilament H (NfH) Antibody
Catalog Number	HM1247
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant NfH protein
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM1250, HM1252 or HM1253 (capture)
Preservatives	0.05% Sodium Azide