

Product Datasheet

Alpha-fetoprotein (AFP) Antibody HM1258

Product Name	Alpha-fetoprotein (AFP) Antibody
Catalog Number	HM1258
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	AFP
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	AFP
Matched Pairs	HM1257
Preservatives	0.1% Sodium Azide