

## Product Datasheet

### Human Epididymis Protein 4 (HE4) Antibody HM1277

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Human Epididymis Protein 4 (HE4) Antibody                                 |
| <b>Catalog Number</b>      | HM1277  |
| <b>Species/Host</b>        | Mouse   |
| <b>Tested Applications</b> | ELISA   |
| <b>Immunogen</b>           | HE4   |
| <b>Brand</b>               | EastCoast Bio   |
| <b>Form/Appearance</b>     | Protein A Purified Monoclonal Antibody                                    |
| <b>Storage</b>             | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| <b>Isotype</b>             | IgG1  |
| <b>Clonality</b>           | Monoclonal  |
| <b>Buffer</b>              | 10 mM Phosphate Buffered Saline, pH 7.2                                   |
| <b>Preservatives</b>       | 0.09% Sodium Azide  |