

Product Datasheet

Human Epididymis Protein 4 (HE4) Antibody HM1278

Product Name	Human Epididymis Protein 4 (HE4) Antibody
Catalog Number	HM1278
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	HE4
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Preservatives	0.09% Sodium Azide