

Product Datasheet

Human Epididymis Protein 4 (HE4) Antibody HM1279

| | |
|----------------------------|---|
| Product Name | Human Epididymis Protein 4 (HE4) Antibody |
| Catalog Number | HM1279 |
| Species/Host | Mouse |
| Tested Applications | ELISA |
| Immunogen | HE4 |
| Brand | EastCoast Bio |
| Form/Appearance | Protein A Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.2 |
| Preservatives | 0.09% Sodium Azide |