

## Product Datasheet

### Thyroglobulin (TG) Antibody HM1264

<b>Product Name</b>	Thyroglobulin (TG) Antibody
<b>Catalog Number</b>	HM1264
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Purified human thyroglobulin
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Specificity</b>	Thyroglobulin
<b>Preservatives</b>	0.09% Sodium Azide