

Product Datasheet

Thyroglobulin (TG) Antibody HM1265

Product Name	Thyroglobulin (TG) Antibody
Catalog Number	HM1265
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Purified human thyroglobulin
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Thyroglobulin
Preservatives	0.09% Sodium Azide