

Product Datasheet

Cortisol, Feline Antibody HM1256

Product Name	Cortisol, Feline Antibody
Catalog Number	HM1256
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Cortisol-KLH
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: ≤-20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl
Specificity	Cortisol
Preservatives	0.05% Sodium Azide