

Product Datasheet

Procalcitonin (PCT) Antibody, Canine HM1205

Product Name	Procalcitonin (PCT) Antibody, Canine
Catalog Number	HM1205
Species/Host	Mouse
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM1203, HM1204 or HM1206
Preservatives	0.09% Sodium Azide