

Product Datasheet

Progesterone, Feline Antibody HM1336

Product Name	Progesterone, Feline Antibody
Catalog Number	HM1336
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Progesterone-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Specificity	Progesterone
Matched Pairs	LA330 (capture)
Preservatives	0.05% Sodium Azide