

Product Datasheet

Testosterone, Feline Antibody HM1302

Product Name	Testosterone, Feline Antibody
Catalog Number	HM1302
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Testosterone-KLH
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: ≤-20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Specificity	Testosterone
Matched Pairs	LA329 (capture)
Preservatives	0.05% Sodium Azide