

Product Datasheet

Dinitrophenyl (DNP) Antibody HM1328

Product Name	Dinitrophenyl (DNP) Antibody
Catalog Number	HM1328
Species/Host	Mouse
Tested Applications	ELISA, Lateral Flow
Immunogen	DNP-KLH
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	PBS, pH 7.0
Preservatives	0.02% Sodium Azide