

## Product Datasheet

### K2-BSA Antigen 170082

<b>Product Name</b>	K2-BSA Antigen
<b>Catalog Number</b>	170082
<b>Description</b>	Synthetic K2, BSA conjugated
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid; Antigen-BSA conjugate
<b>Concentration</b>	Lot-specific
<b>Storage</b>	Short term storage ( $\leq 2$ weeks): 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
<b>Note</b>	This product was previously available with a Catalog code P83-99-51A and LA288. It has been produced using the same conjugation protocol.
<b>Buffer</b>	Distilled water
<b>Matched Pairs</b>	JP098, JP099 or JP100 (detection)
<b>Shelf Life</b>	60 months