

Product Datasheet

Anti-h BNP 11904 SPTN-5 100970

Product Name	Anti-h BNP 11904 SPTN-5
Catalog Number	100970
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Alternative Names	B-type natriuretic peptide
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	C-terminal epitope in amino acid region 26-32.
Purity	≥ 95 %
Associated Products	Synthetic Brain Natriuretic Peptide-32 (BNP), Medix Biochemica 129-10 and synthetic BNP antigen, Medix Biochemica LA788.
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.1-6.6
Cross Reactivity	Does not recognize ANP or CNP.
Specificity	Antibody recognizes human B-type natriuretic peptide
Matched Pairs	BNP 11906 and 11908
Shelf Life	Unspecified