

Product Datasheet

Anti-h BNP 11906 SPTN-5 100971

Product Name	Anti-h BNP 11906 SPTN-5
Catalog Number	100971
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Alternative Names	B-type natriuretic peptide
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D (N-terminal or mid-region based on pairing information)
Purity	≥ 95 %
Associated Products	Synthetic Brain Natriuretic Peptide-32 (BNP), Medix Biochemica 129-10 and synthetic BNP antigen, Medix Biochemica LA788.
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	7.2-8.0
Cross Reactivity	Does not recognize ANP or CNP
Specificity	Antibody recognizes human B-type natriuretic peptide
Matched Pairs	BNP 11904
Shelf Life	Unspecified