

Product Datasheet

D-dimer Antibody HM1408

Product Name	D-dimer Antibody
Catalog Number	HM1408
Species/Host	in vitro cell culture
Tested Applications	ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG1
Clonality	Monoclonal
Buffer	PBS, pH 7.4
Specificity	D-dimer, D-dimer β -chain, and fibrin degradation products
Matched Pairs	HM1409 (detection)
Preservatives	0.09% Sodium Azide