

Product Datasheet

Donkey anti-Sheep IgG (H&L) Antibody JP203

Product Name	Donkey anti-Sheep IgG (H&L) Antibody
Catalog Number	JP203
Species/Host	Donkey
Tested Applications	ELISA
Immunogen	Purified Sheep IgG obtained from normal sheep serum
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Polyclonal Antibody
Storage	2-8°C
Isotype	IgG
Clonality	Polyclonal
Purity	>95% by SDS-PAGE
Buffer	PBS, pH 7.2
Specificity	H & L chains Sheep IgG
Preservatives	0.05% Sodium Azide