

Product Datasheet

Ethyl glucuronide (EtG) Antibody HM1400

Product Name	Ethyl glucuronide (EtG) Antibody
Catalog Number	HM1400
Species/Host	Sheep
Tested Applications	ELISA, Lateral Flow
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Purity	>95%
Buffer	PBS, pH 7.4
Specificity	Ethyl glucuronide
Matched Pairs	LA809 (capture)