

Product Datasheet

Ethyl glucuronide (EtG)-BSA Antigen LA809

Product Name	Ethyl glucuronide (EtG)-BSA Antigen
Catalog Number	LA809
Tested Applications	ELISA, Lateral Flow
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	PBS, pH 7.4
Matched Pairs	HM1400 (detection)