

Product Datasheet

Monkeypox Virus (MPXV-A29L) Antibody HM1338

Product Name	Monkeypox Virus (MPXV-A29L) Antibody
Catalog Number	HM1338
Species/Host	in vitro cell culture
Tested Applications	ELISA, Lateral Flow
Immunogen	Purified A29L protein
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	Tris Glycine, pH 7.0-7.4
Matched Pairs	HM1339 (detection)