

## Product Datasheet

### Monkeypox Virus (MPXV-A29L) Antigen, Recombinant >95% LA808

**Product Name** Monkeypox Virus (MPXV-A29L) Antigen, Recombinant >95%

**Catalog Number** LA808

**Tested Applications** ELISA, Lateral Flow

**Brand** EastCoast Bio

**Form/Appearance** Affinity Purified Recombinant Antigen (His-tag)

**Storage** Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Epitope** A29L

**Purity** >95% by SDS-PAGE

**Buffer** 50 mM Tris, 100 mM Glycine, pH 7.0

**Matched Pairs** HM1338, HM1339

**Expression System** Mammalian Cell Expression System (HEK293)