

Product Datasheet

Ubiquitin carboxy terminal Hydrolase (UCHL1) Antibody HM1418

Product Name	Ubiquitin carboxy terminal Hydrolase (UCHL1) Antibody
Catalog Number	HM1418
Species/Host	Mouse
Tested Applications	ELISA, Lateral Flow
Immunogen	Human UCLH1 (HEK293 expressed)
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM1420, HM1421 or HM1422 (detection)
Preservatives	0.05% Sodium Azide