

Product Datasheet

Anti-Human CD62E Monoclonal Antibody, Azide Free Clone B-P7 852.730

Product Name	Anti-Human CD62E Monoclonal Antibody, Azide Free Clone B-P7
Catalog Number	852.730
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Activated human endothelial cells
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Lambda light chain
Clonality	Monoclonal
Specificity	Recognises the Endothelial Leucocyte Adhesion Molecule-1 (ELAM-1), E-selectin, a 97 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells