

## Product Datasheet

### Anti-Human CD102 Monoclonal Antibody, Unconjugated Clone B-T1 852.763.020

<b>Product Name</b>	Anti-Human CD102 Monoclonal Antibody, Unconjugated Clone B-T1
<b>Catalog Number</b>	852.763.020
<b>Description</b>	Mouse monoclonal
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Purified
<b>Tested Applications</b>	IP, FC
<b>Immunogen</b>	ICAM-2 transfected CHO cell line
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
<b>Storage</b>	Stable at +2-8°C for 6 months. For longer storage freeze aliquots.
<b>Note</b>	Use 10 µl to label 10 <sup>6</sup> cells or 100 µl of whole blood
<b>Isotype</b>	Mouse IgG1 Kappa light chain
<b>Clonality</b>	Monoclonal
<b>Specificity</b>	Recognises the Intercellular Adhesion Molecule-2 (ICAM-2), a 55-65 kDa protein
<b>Hybridoma</b>	Myeloma X63/AG.8653 x Balb/c spleen cells